

产品名称: C7 Rabbit Polyclonal Antibody
产品货号: APRab01369



产品概述 (Summary)

产品名称 (Production Name)	C7 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	C7
别名 (Alternative Names)	Complement component C7
基因 ID (Gene ID)	730
蛋白 ID (SwissProt ID)	P10643.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 94 kDa; Observed MW: 94 kDa

研究背景 (Background)

Constituent of the membrane attack complex (MAC) that plays a key role in the innate and adaptive immune response by

产品名称: C7 Rabbit Polyclonal Antibody
产品货号: APRab01369




forming pores in the plasma membrane of target cells.


研究领域 (Research Area)

Immunology


图片 (Image Data)



Western blot analysis of C7 in C6, PC-12, 3T3, Raw264.7 lysates using C7 antibody.



Western blot analysis of C7 in HeLa, 293 lysates using C7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using C7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: C7 Rabbit Polyclonal Antibody
产品货号: APRab01369



注意事项 (Note)

For research use only .