

产品名称: CAP1 Rabbit Polyclonal Antibody
产品货号: APRab01371



产品概述 (Summary)

产品名称 (Production Name)	CAP1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CAP1
别名 (Alternative Names)	CAP; CAP1-PEN
基因 ID (Gene ID)	10487
蛋白 ID (SwissProt ID)	Q01518.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 52 kDa; Observed MW: 56 kDa

研究背景 (Background)

Directly regulates filament dynamics and has been implicated in a number of complex developmental and morphological

产品名称: CAP1 Rabbit Polyclonal Antibody
产品货号: APRab01371




processes, including mRNA localization and the establishment of cell polarity.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of CAP1 in 293T, MCF-7 lysates using CAP1 antibody.



Immunocytochemistry analysis of CAP1(green) in HeLa using CAP1 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using CAP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: CAP1 Rabbit Polyclonal Antibody
产品货号: APRab01371



注意事项 (Note)

For research use only .