

产品名称: SENP1 Rabbit Polyclonal Antibody
产品货号: APRab01373



产品概述 (Summary)

产品名称 (Production Name)	SENP1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SENP1
别名 (Alternative Names)	SENP1; SuPr2
基因 ID (Gene ID)	29843
蛋白 ID (SwissProt ID)	Q9P0U3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 73 kDa; Observed MW: 76 kDa

研究背景 (Background)

Protease that catalyzes two essential functions in the SUMO pathway: processing of full-length SUMO1, SUMO2 and

产品名称: SENP1 Rabbit Polyclonal Antibody
产品货号: APRab01373




SUMO3 to their mature forms and deconjugation of SUMO1, SUMO2 and SUMO3 from targeted proteins.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Western blot analysis of SENP1 in 293, HeLa, COS7, PC-3 lysates using SENP1 antibody.



Immunocytochemistry analysis of SENP1(green) in HeLa using SENP1 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using SENP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: SENP1 Rabbit Polyclonal Antibody
产品货号: APRab01373



注意事项 (Note)

For research use only .