

产品名称: Ubiquitin 60S Ribosomal Protein L40 Rabbit Polyclonal Antibody
产品货号: APRab01374



产品概述 (Summary)

产品名称 (Production Name)	Ubiquitin 60S Ribosomal Protein L40 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC, IP
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UBA52
别名 (Alternative Names)	L40; CEP52; RPL40; HUBCEP52
基因 ID (Gene ID)	7311
蛋白 ID (SwissProt ID)	P62987.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, FC 1:50-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 15 kDa; Observed MW: 10 kDa

研究背景 (Background)

产品名称: Ubiquitin 60S Ribosomal Protein L40 Rabbit
Polyclonal Antibody
产品货号: APRab01374




Ubiquitin: Exists either covalently attached to another protein, or free (unanchored).


研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Western blot analysis of UBA52 in C6 lysates using UBA52 antibody.



Western blot analysis of UBA52 in HeLa, SKOV-3 lysates using UBA52 antibody.



Immunocytochemistry analysis of UBA52(green) in HeLa using UBA52 antibody, and DAPI(blue)

**产品名称: Ubiquitin 60S Ribosomal Protein L40 Rabbit
Polyclonal Antibody**
产品货号: APRab01374



注意事项 (Note)

For research use only .