

产品名称: GRK3 Rabbit Polyclonal Antibody
产品货号: APRab01375



产品概述 (Summary)

产品名称 (Production Name)	GRK3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRK3
别名 (Alternative Names)	ADRBK2; BARK2; GRK3
基因 ID (Gene ID)	157
蛋白 ID (SwissProt ID)	P35626.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 80 kDa; Observed MW: 80 kDa

研究背景 (Background)

Specifically phosphorylates the agonist-occupied form of the beta-adrenergic and closely related receptors.


产品名称: GRK3 Rabbit Polyclonal Antibody
产品货号: APRab01375




研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of GRK3 in Jurkat lysates using GRK3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using GRK3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .