

产品名称: Integrin beta 1 Rabbit Polyclonal Antibody
产品货号: APRab01383



产品概述 (Summary)

产品名称 (Production Name)	Integrin beta 1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ITGB1
别名 (Alternative Names)	ITGB1; FNRB; MDF2; MSK12; Integrin beta-1; Fibronectin receptor subunit beta; VLA-4 subunit beta; CD antigen CD29
基因 ID (Gene ID)	3688
蛋白 ID (SwissProt ID)	P05556.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 88 kDa; Observed MW: 120-160 kDa

研究背景 (Background)

产品名称: Integrin beta 1 Rabbit Polyclonal Antibody
产品货号: APRab01383




Integrins are α/β heterodimeric cell surface receptors that play a pivotal role in cell adhesion and migration, as well as in growth and survival. The integrin family contains at least 18 α and 8 β subunits that form 24 known integrins with distinct tissue distribution and overlapping ligand specificities.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Integrin beta 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .