

产品名称: GRB2 Rabbit Polyclonal Antibody
产品货号: AP Rab01386



产品概述 (Summary)

产品名称 (Production Name)	GRB2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRB2
别名 (Alternative Names)	GRB2; ASH; Growth factor receptor-bound protein 2; Adapter protein GRB2; Protein Ash; SH2/SH3 adapter GRB2
基因 ID (Gene ID)	2885
蛋白 ID (SwissProt ID)	P62993.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 25 kDa; Observed MW: 25 kDa

研究背景 (Background)

产品名称: GRB2 Rabbit Polyclonal Antibody
产品货号: APRab01386




Adapter protein that provides a critical link between cell surface growth factor receptors and the Ras signaling pathway.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of GRB2 in PC-12, C6 lysates using GRB2 antibody.



Western blot analysis of GRB2 in MCF-7, 3T3, Raw264.7, HeLa lysates using GRB2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using GRB2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: GRB2 Rabbit Polyclonal Antibody
产品货号: APRab01386



注意事项 (Note)

For research use only .