

产品名称: Phospho-p38 (Thr180/Tyr182) Rabbit Polyclonal Antibody
产品货号: APRab01390



产品概述 (Summary)

产品名称 (Production Name)	Phospho-p38 (Thr180/Tyr182) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAPK14 MAPK14; CSBP; CSBP1; CSBP2; CSPB1; MXI2; SAPK2A; Mitogen-activated protein kinase 14; MAP kinase 14; MAPK 14; Cytokine suppressive anti-inflammatory drug-binding protein; CSAID-binding protein; CSBP; MAP kinase MXI2; MAX-interacting protein
别名 (Alternative Names)	
基因 ID (Gene ID)	1432
蛋白 ID (SwissProt ID)	Q16539.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
----------------------	-----------------

产品名称: Phospho-p38 (Thr180/Tyr182) Rabbit Polyclonal Antibody
产品货号: APRab01390



蛋白分子量 (Molecular Weight) Calculated MW: 41 kDa; Observed MW: 41 kDa


研究背景 (Background)

Responds to activation by environmental stress, pro-inflammatory cytokines and lipopolysaccharide (LPS) by phosphorylating a number of transcription factors, such as ELK1 and ATF2 and several downstream kinases, such as MAPKAPK2 and MAPKAPK5. Plays a critical role in the production of some cytokines, for example IL-6.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of p38 (Phospho- T180/Y182) in 293, 3T3 lysates using Phospho-p38 (Thr180/Tyr182) antibody.

注意事项 (Note)

For research use only .