

产品名称: Fascin Rabbit Polyclonal Antibody
产品货号: APRab01395



产品概述 (Summary)

产品名称 (Production Name)	Fascin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FSCN1 55 kDa actin bundling protein; Actin bundling protein; FAN1; Fascin 1; Fascin;
别名 (Alternative Names)	Singed (Drosophila) like (sea urchin fascin homolog like); Fascin homolog 1; Fascin
基因 ID (Gene ID)	6624
蛋白 ID (SwissProt ID)	Q16658.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 55 kDa; Observed MW: 55 kDa

产品名称: Fascin Rabbit Polyclonal Antibody
产品货号: APRab01395




研究背景 (Background)

Promotes cross-linkage of parallel actin filaments during the formation of cell protrusions (lamellipodia and filopodia), and therefore plays an important role in regulating cell migration. It has been reported that fascin may also regulate filopodia formation by a mechanism independent of its actin-bundling functions, though less is known about this mechanism of action.


研究领域 (Research Area)

Cell Biology

图片 (Image Data)




Immunocytochemistry analysis of Fascin (green) in 3T3 using Fascin antibody, and DAPI (blue).




Western blot analysis of Fascin in mouse skin, rat skin lysates using Fascin antibody.

产品名称: Fascin Rabbit Polyclonal Antibody
产品货号: APRab01395



Western blot analysis of Fascin in HeLa, A549, HepG2, COS7 lysates using Fascin antibody



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Fascin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .