

产品名称: Platin L Rabbit Polyclonal Antibody
产品货号: APRab01400



产品概述 (Summary)

产品名称 (Production Name)	Platin L Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LCP1
别名 (Alternative Names)	CP64; L platin; LC64P; LCP1; LPL; Lplatin; Platin 2; PLS2
基因 ID (Gene ID)	3936
蛋白 ID (SwissProt ID)	P13796.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 70 kDa; Observed MW: 70 kDa

研究背景 (Background)

Actin-binding protein. Plays a role in the activation of T-cells in response to costimulation through TCR/CD3 and CD2 or

产品名称: Platin L Rabbit Polyclonal Antibody
产品货号: AP Rab01400




CD28. Modulates the cell surface expression of IL2RA/CD25 and CD69.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Immunocytochemistry analysis of Platin L (green) in Jurkat using Platin L antibody, and DAPI (blue).



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Platin L antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .