

产品名称: Human IgG Rabbit Polyclonal Antibody
产品货号: AP Rab03306



产品概述 (Summary)

产品名称 (Production Name)	Human IgG Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IGHG1
别名 (Alternative Names)	IGHG1; Ig gamma-1 chain C region
基因 ID (Gene ID)	-
蛋白 ID (SwissProt ID)	-.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	-

研究背景 (Background)

IgG is a monomeric immunoglobulin, built of two heavy chains gamma and two light chains. Each molecule has two

产品名称: Human IgG Rabbit Polyclonal Antibody
产品货号: APRab03306




antigen binding sites. This is the most abundant immunoglobulin and is approximately equally distributed in blood and in tissue liquids, constituting 75% of serum immunoglobulins in humans.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Indirect ELISA assay for antiHuman IgG rabbit pAb. Antigen coating concentration: 4ug/ml.

注意事项 (Note)

For research use only .