

产品名称: NFAT2 Rabbit Polyclonal Antibody
产品货号: APRab03352



产品概述 (Summary)

产品名称 (Production Name)	NFAT2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NFATC1 NFATC1; NFAT2; NFATC; Nuclear factor of activated T-cells; cytoplasmic 1; NF-
别名 (Alternative Names)	ATc1; NFATc1; NFAT transcription complex cytosolic component; NF-ATc; NFATc
基因 ID (Gene ID)	4772
蛋白 ID (SwissProt ID)	O95644.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 101 kDa; Observed MW: 105 kDa

产品名称: NFAT2 Rabbit Polyclonal Antibody
产品货号: APRab03352




研究背景 (Background)

NFAT proteins are predominantly expressed in cells of the immune system, but are also expressed in skeletal muscle, keratinocytes, and adipocytes, regulating cell differentiation programs in these cells.


研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of NFAT2 in HuvEc lysates using NFAT2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung using NFAT2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .