

**产品名称:** DRP1 Rabbit Polyclonal Antibody

**产品货号:** APRab03371

### 产品概述 (Summary)

|                   |                                                                                                |
|-------------------|------------------------------------------------------------------------------------------------|
| 描述(Description)   | Rabbit polyclonal Antibody                                                                     |
| 宿主(Host)          | Rabbit                                                                                         |
| 应用(Application)   | WB,IHC                                                                                         |
| 种属反应(Reactivity)  | Human,Mouse,Rat                                                                                |
| 偶联物(Conjugation)  | Unconjugated                                                                                   |
| 修饰(Modification)  | Unmodified                                                                                     |
| 同种型(Isotype)      | IgG                                                                                            |
| 克隆(Clonality)     | Polyclonal                                                                                     |
| 剂型(Form)          | Liquid                                                                                         |
| 保存条件(Storage)     | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.       |
| 储存溶液(Buffer)      | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3. |
| 纯化 (Purification) | Affinity Purification                                                                          |

### 产品应用(Application)

|                       |                                            |
|-----------------------|--------------------------------------------|
| 稀释比(Dilution Ratio)   | WB 1:500-1:1000,IHC 1:50-1:100             |
| 分子量(Molecular Weight) | Calculated MW: 82 kDa; Observed MW: 82 kDa |

### 抗原信息(Antigen Information)

|                       |                                                                                                                                                                                                                                                  |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名(Gene Name)        | DNM1L                                                                                                                                                                                                                                            |
| 别名(Alternative Names) | DNM1L; DLP1; DRP1; Dynamin-1-like protein; Dnm1p/Vps1p-like protein; DVLP; Dynamin family member proline-rich carboxyl-terminal domain less; Dymple; Dynamin-like protein; Dynamin-like protein 4; Dynamin-like protein IV; HdynIV; Dynamin-rela |
| 基因 ID(Gene ID)        | 10059                                                                                                                                                                                                                                            |
| SwissProt ID          | O00429                                                                                                                                                                                                                                           |
| 免疫原(Immunogen)        |                                                                                                                                                                                                                                                  |

### 研究背景 (Background)

The protein encoded by this gene is a member of the dynamin superfamily of GTPases. Members of the dynamin-related

subfamily, including the *S. cerevisiae* proteins Dnm1 and Vps1, contain the N-terminal tripartite GTPase domain but do not have the pleckstrin homology or proline-rich domains. This protein establishes mitochondrial morphology through a role in distributing mitochondrial tubules throughout the cytoplasm. The gene has 3 alternatively spliced transcripts encoding different isoforms. These transcripts are alternatively polyadenylated.


## 研究领域 (Research Area)

Neuroscience


## 注意事项 (Note)

For research use only.


## 图片 (Image Data)




Western blot analysis of DRP1 in HeLa lysates using DRP1 antibody.



Immunofluorescence analysis of DRP1 in rat lung using DRP1 antibody (red), and DAPI (blue).



Immunofluorescence analysis of DRP1 in mouse kidney using DRP1 antibody (red), and DAPI (blue).



Immunohistochemistry analysis of paraffin-embedded Human stomach cancer using DRP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.

## EnkiLife 优势产品推荐

| 产品类别                        | 产品货号     | 产品名称                                               |
|-----------------------------|----------|----------------------------------------------------|
| WB 解决方案                     | RA10020  | 2h 极速 WB 即用型全流程试剂盒                                 |
|                             | RA10021  | 4h 快速 WB 即用型全流程试剂盒                                 |
|                             | RA10042  | 5h 畅享版 WB 全流程试剂盒                                   |
|                             | RA10037  | 校准级彩色预染蛋白 Marker (8-180kDa)                        |
|                             | RA10038  | 校准级彩色预染蛋白 Marker (10-250kDa)                       |
|                             | RA10039  | 校准级高分子彩色预染蛋白 Marker (25-400kDa)                    |
| TSA 多重荧光染色试剂盒               | RA10008  | TSA 双标三色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10009  | TSA 三标四色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10010  | TSA 四标五色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10011  | TSA 五标六色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10012  | TSA 六标七色多重荧光染色试剂盒 (mIHC)                           |
| IHC 检测试剂盒                   | RA10006  | HRP Anti-Mouse/Rabbit IHC Detection System         |
|                             | RA10007  | Polymer-HRP Anti-Mouse/Rabbit IHC Detection System |
| 抗体标记试剂盒                     | RE80004p | 辣根过氧化物酶(HRP)抗体标记试剂盒                                |
|                             | RE80002q | Sulfo-NHS-生物素标记试剂盒                                 |
|                             | RE80007p | Cy3 荧光素标记试剂盒                                       |
|                             | RE80011p | Fluor488 荧光素标记试剂盒                                  |
|                             | RE80017p | Fluor750 荧光素标记试剂盒                                  |
|                             | RE80005p | 藻红蛋白(R-PE) 抗体快速标记试剂盒                               |
|                             | RE80040  | PE-Cy7 串联染料抗体快速标记试剂盒                               |
| 稳转细胞系构建服务<br>(免费赠送全膜 WB 验证) | TS-0001  | 过表达稳转细胞系构建                                         |
|                             | TS-0002  | 敲低稳转细胞系构建                                          |
|                             | TS-0003  | 敲除细胞系构建                                            |

限于篇幅，更多产品详见官网：<https://www.enkilife.cn>