

产品名称: PPAR alpha Rabbit Polyclonal Antibody
产品货号: AP Rab03383



产品概述 (Summary)

产品名称 (Production Name)	PPAR alpha Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPARA
别名 (Alternative Names)	PPARA; NR1C1; PPAR; Peroxisome proliferator-activated receptor alpha; PPAR-alpha; Nuclear receptor subfamily 1 group C member 1
基因 ID (Gene ID)	5465
蛋白 ID (SwissProt ID)	Q07869.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 52 kDa; Observed MW: 52 kDa

研究背景 (Background)

产品名称: PPAR alpha Rabbit Polyclonal Antibody
产品货号: APRab03383




Ligand-activated transcription factor.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Western blot analysis of PPAR alpha in 3T3 lysates using PPAR alpha antibody.



Immunohistochemistry analysis of paraffin-embedded mouse spleen tissue using PPARα antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of PPAR alpha in HeLa cell using PPARα antibody (red) and Galectin3 antibody (green)

产品名称: PPAR alpha Rabbit Polyclonal Antibody
产品货号: APRab03383



注意事项 (Note)

For research use only .