

产品名称: RPL7 Rabbit Polyclonal Antibody
产品货号: APRab03399



产品概述 (Summary)

产品名称 (Production Name)	RPL7 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPL7
别名 (Alternative Names)	L7; humL7-1
基因 ID (Gene ID)	6129
蛋白 ID (SwissProt ID)	P18124.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 29 kDa; Observed MW: 32 kDa

研究背景 (Background)

Component of the large ribosomal subunit (PubMed:12962325). Binds to G-rich structures in 28S rRNA and in mRNAs.

产品名称: RPL7 Rabbit Polyclonal Antibody
产品货号: APRab03399




Plays a regulatory role in the translation apparatus; inhibits cell-free translation of mRNAs.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Western blot analysis of RPL7 in various lysates using RPL7 antibody.




Western blot analysis of RPL7 in cat sample lysates using Ribosomal Protein L7 antibody.




Western blot analysis of RPL7 in mouse kidney, mouse brain, HELA, KB, SH-SY5Y, 293T and 3T3 lysates using Ribosomal Protein L7 antibody.

产品名称: RPL7 Rabbit Polyclonal Antibody
产品货号: APRab03399



Western blot analysis of RPL7 in 293 and HeLa lysates using RPL7 antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of PL7 in K562 using RPL7 antibody.

注意事项 (Note)

For research use only .