

产品名称: 14-3-3 Rabbit Polyclonal Antibody
产品货号: APRab03629



产品概述 (Summary)

产品名称 (Production Name)	14-3-3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	YWHAB/YWHAЕ/YWHAH 14-3-3; 14 3 3; 14 3 3 protein beta; 14 3 3 protein beta/alpha; 14 3 3 protein zeta; 14 3 3 zeta; 14-3-3 protein beta/alpha; 14-3-3 protein/cytosolic phospholipase A2; 1433B_HUMAN; GW128; HS1; KCIP 1; KCIP-1; MGC111427; MGC126532; MGC138156; N-terminally processed; Protein 1054; Protein kinase C inhibitor protein 1; Tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein; delta polypeptide; Tyrosine 3/tryptophan 5 -monooxygenase activation protein; zeta polypeptide; YWHAB; YWHAD; YWHAZ.
别名 (Alternative Names)	
基因 ID (Gene ID)	7529/7531/7533/7532/2810/10971/7534
蛋白 ID (SwissProt ID)	P31946/P62258/Q04917/P61981/P31947/P27348/P63104.

产品名称: 14-3-3 Rabbit Polyclonal Antibody
产品货号: APRab03629



产品应用 (Application)

稀释比 (Dilution Ratio)

WB 1:500-1:1000

蛋白分子量 (Molecular Weight)

Calculated MW: 27-29 kDa; Observed MW: 27-29 kDa


研究背景 (Background)

Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathways. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner. Negative regulator of osteogenesis.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of 14-3-3 (pan) in rat Hippocampus, rat Cerebral Cortex and rat Brain lysates using 14-3-3 (pan) antibody.



Western blot analysis of 14-3-3 in C6, 3T3 and Hela lysates using 14-3-3 (pan) antibody.

产品名称: 14-3-3 Rabbit Polyclonal Antibody
产品货号: APRab03629



注意事项 (Note)

For research use only .