

产品名称: GABARAP Rabbit Polyclonal Antibody
产品货号: APRab03631



产品概述 (Summary)

产品名称 (Production Name)	GABARAP Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GABARAP
别名 (Alternative Names)	GABARAP; GABARAPL1; GABARAPL2
基因 ID (Gene ID)	11337
蛋白 ID (SwissProt ID)	O95166.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 14 kDa; Observed MW: 14, 16 kDa

研究背景 (Background)

May play a role in intracellular transport of GABA(A) receptors and its interaction with the cytoskeleton.


产品名称: GABARAP Rabbit Polyclonal Antibody
产品货号: APRab03631




研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of GABARAP in A549, 3T3 and HeLa lysates using GABARAP antibody.



Western blot analysis of GABARAP in rat Hippocampus, rat Cerebral Cortex, rat Brain and A549 lysates using GABARAP antibody.

注意事项 (Note)

For research use only .