

产品名称: FOXP3 Rabbit Polyclonal Antibody
产品货号: APRab03820



产品概述 (Summary)

产品名称 (Production Name)	FOXP3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FOXP3
别名 (Alternative Names)	FOXP3; IPEX; JM2; Forkhead box protein P3; Scurfin
基因 ID (Gene ID)	50943
蛋白 ID (SwissProt ID)	Q9BZS1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 47 kDa

研究背景 (Background)

Defects in FOXP3 are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX)

产品名称: FOXP3 Rabbit Polyclonal Antibody
产品货号: APRab03820




[MIM:304790]; also known as X-linked autoimmunity-immunodeficiency syndrome. IPEX is characterized by neonatal onset insulin-dependent diabetes mellitus, infections, secretory diarrhea, trombocytopenia, anemia and eczema. It is usually lethal in infancy.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Western blot analysis of FOXP3 in PC-12 lysates using FOXP3 antibody.




Immunohistochemistry analysis of paraffin-embedded Human tonsilla using FOXP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.




产品名称: FOXP3 Rabbit Polyclonal Antibody
产品货号: APRab03820




Immunohistochemistry analysis of paraffin-embedded Human bladder cancer using FOXP3 antibody. High-pressure and temperature Tris-EDTA was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human colon using FOXP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human colon using FOXP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human colon using FOXP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

产品名称: FOXP3 Rabbit Polyclonal Antibody
产品货号: APRab03820



For research use only .