

产品名称: Smad1/5/9 Rabbit Polyclonal Antibody
产品货号: AP Rab03822



产品概述 (Summary)

产品名称 (Production Name)	Smad1/5/9 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SMAD1/SMAD5/SMAD9 Mothers against decapentaplegic homolog 1/5/9; MAD homolog 1; Mothers against DPP homolog 1; JV4-1; Mad-related protein 1; SMAD family member 1; SMAD 1/5/9; Smad1; hSMAD1; Transforming growth factor-beta-signaling protein 1; BSP-1
别名 (Alternative Names)	
基因 ID (Gene ID)	4086/4090/4093
蛋白 ID (SwissProt ID)	Q15797/Q99717/O15198.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 52 kDa; Observed MW: 52 kDa

产品名称: Smad1/5/9 Rabbit Polyclonal Antibody
产品货号: APRab03822




研究背景 (Background)

Transcriptional modulator activated by BMP (bone morphogenetic proteins) type 1 receptor kinase. SMAD1 is a receptor-regulated SMAD (R-SMAD). SMAD1/OAZ1/PSMB4 complex mediates the degradation of the CREBBP/EP300 repressor SNIP1. May act synergistically with SMAD4 and YY1 in bone morphogenetic protein (BMP)-mediated cardiac-specific gene expression.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Smad1/5/9 in 293T lysates using Smad1/5/9 antibody.

注意事项 (Note)

For research use only .