

产品名称: Cyclin E1 Rabbit Polyclonal Antibody
产品货号: APRab03827



产品概述 (Summary)

产品名称 (Production Name)	Cyclin E1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CCNE1
别名 (Alternative Names)	CCNE1; CCNE; G1/S-specific cyclin-E1
基因 ID (Gene ID)	898
蛋白 ID (SwissProt ID)	P24864.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 47 kDa

研究背景 (Background)

Essential for the control of the cell cycle at the G1/S (start) transition. a member of the highly conserved cyclin family,

产品名称: Cyclin E1 Rabbit Polyclonal Antibody
产品货号: APRab03827




whose members are characterized by a dramatic periodicity in protein abundance through the cell cycle. Cyclins function as regulators of CDK kinases.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Cyclin E1 in 3T3 and HeLa lysates using Cyclin E1 antibody.

注意事项 (Note)

For research use only .