

产品名称: ERK1/2 Rabbit Polyclonal Antibody
产品货号: APRab03846



产品概述 (Summary)

产品名称 (Production Name)	ERK1/2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAPK3/MAPK1 MAPK3; ERK1; ERT2; ERK-1; PRKM3; P44ERK1; P44MAPK; HS44KDAP;
别名 (Alternative Names)	HUMKER1A; p44-ERK1; p44-MAPK; MAPK1; ERK; p38; p40; p41; ERK2; ERT1; ERK-2; MAPK2; PRKM1; PRKM2; P42MAPK; p41mapk; p42-MAPK.
基因 ID (Gene ID)	5595/5594
蛋白 ID (SwissProt ID)	P27361/P28482.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 42,44 kDa; Observed MW: 42,44 kDa

产品名称: ERK1/2 Rabbit Polyclonal Antibody
产品货号: AP Rab03846




研究背景 (Background)

Serine/threonine kinase which acts as an essential component of the MAP kinase signal transduction pathway. MAPK1/ERK2 and MAPK3/ERK1 are the 2 MAPKs which play an important role in the MAPK/ERK cascade. They participate also in a signaling cascade initiated by activated KIT and KITLG/SCF. Depending on the cellular context, the MAPK/ERK cascade mediates diverse biological functions such as cell growth, adhesion, survival and differentiation through the regulation of transcription, translation, cytoskeletal rearrangements.


研究领域 (Research Area)

Cell Biology

图片 (Image Data)




Western blot analysis of ERK1/2 in COLO205 lysates using ERK1/2 antibody.



Western blot analysis of ERK1/2 in K562 lysates using ERK 1/2 antibody.

产品名称: ERK1/2 Rabbit Polyclonal Antibody
产品货号: APRab03846



Western blot analysis of ERK1/2 in HELA lysates using ERK 1/2 antibody.

注意事项 (Note)

For research use only .