

产品名称: Phospho-Histone H3 (Ser10) Rabbit Polyclonal Antibody
产品货号: APRab03869



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Histone H3 (Ser10) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	H3C1
别名 (Alternative Names)	H3 histone; family 3A; H3 histone; family 3B (H3.3B); H3.3A; H3.3B; H33; H3F3; H3F3A; H3F3B; Histone H3.3
基因 ID (Gene ID)	8350
蛋白 ID (SwissProt ID)	P68431.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 15 kDa; Observed MW: 15 kDa

产品名称: Phospho-Histone H3 (Ser10) Rabbit Polyclonal Antibody
产品货号: APRab03869




研究背景 (Background)

H3 Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of Phospho-Histone H3 (Ser10) in HeLa lysates, untreated or 1mM hydroxyurea treated for 5 hours, using Phospho-Histone H3(Ser10) antibody.



Immunofluorescence analysis of Phospho-Histone H3 (Ser10) in HeLa cells labeled using Phospho-Histone H3(Ser10) antibody, and DAPI (blue).

注意事项 (Note)

For research use only .