

产品名称: 14-3-3 β/ζ (phospho Ser184/186) Rabbit Polyclonal Antibody
产品货号: APRab04176



产品概述 (Summary)

产品名称 (Production Name)	14-3-3 β/ζ (phospho Ser184/186) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	YWHAB/YWHAZ
别名 (Alternative Names)	YWHAZ; 14-3-3 protein zeta/delta; Protein kinase C inhibitor protein 1; KCIP-1; YWHAB; 14-3-3 protein beta/alpha; Protein 1054; Protein kinase C inhibitor protein 1; KCIP-1
基因 ID (Gene ID)	7534/7529
蛋白 ID (SwissProt ID)	P63104/P31946.The antiserum was produced against synthesized peptide derived from human 14-3-3 beta/zeta around the phosphorylation site of Ser184/186. AA range:151-200

产品应用 (Application)

产品名称: 14-3-3 β/ζ (phospho Ser184/186) Rabbit
Polyclonal Antibody
产品货号: APRab04176



稀释比 (Dilution Ratio)


WB 1:500-1:2000, ELISA 1:5000-1:10000

蛋白分子量 (Molecular Weight)

研究背景 (Background)

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from RAW264.7 cells treated with UV 15', using 14-3-3 β/ζ (Phospho-Ser184/186) Antibody. The lane on the right is blocked with the phospho peptide.

注意事项 (Note)

For research use only .