

产品名称: Bcl-2 (phospho Ser70) Rabbit Polyclonal Antibody
产品货号: AP Rab04304



产品概述 (Summary)

产品名称 (Production Name)	Bcl-2 (phospho Ser70) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BCL2
别名 (Alternative Names)	BCL2; Apoptosis regulator Bcl-2
基因 ID (Gene ID)	24224.0
蛋白 ID (SwissProt ID)	P49950.The antiserum was produced against synthesized peptide derived from human BCL2 around the phosphorylation site of Ser70. AA range:36-85

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	26kDa

研究背景 (Background)

产品名称: Bcl-2 (phospho Ser70) Rabbit Polyclonal Antibody
产品货号: APRab04304




BCL2 encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Alternative splicing results in multiple transcript variants.


研究领域 (Research Area)

Stem cell pathway; PI3K/Akt; NF_kappaB


图片 (Image Data)



Western blot analysis of BCL2 (Phospho-Ser70) Antibody. The lane on the right is blocked with the BCL2 (Phospho-Ser70) peptide.



Western Blot analysis of JK cells using Phospho-Bcl-2 (S70) Polyclonal Antibody



产品名称: Bcl-2 (phospho Ser70) Rabbit Polyclonal Antibody
产品货号: APRab04304



Western Blot analysis of Jurkat cells using Phospho-Bcl-2 (S70) Polyclonal Antibody

注意事项 (Note)

For research use only .