

产品名称: GATA-3 (phospho Ser308) Rabbit Polyclonal Antibody
产品货号: APRab04726



产品概述 (Summary)

产品名称 (Production Name)	GATA-3 (phospho Ser308) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GATA3
别名 (Alternative Names)	GATA3; Trans-acting T-cell-specific transcription factor GATA-3; GATA-binding factor 3
基因 ID (Gene ID)	2625.0
蛋白 ID (SwissProt ID)	P23771. The antiserum was produced against synthesized peptide derived from human GATA3 around the phosphorylation site of Ser308. AA range: 274-323

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:20000
----------------------	--

产品名称: GATA-3 (phospho Ser308) Rabbit Polyclonal Antibody
产品货号: APRab04726



蛋白分子量 (Molecular Weight)


研究背景 (Background)

This gene encodes a protein which belongs to the GATA family of transcription factors. The protein contains two GATA-type zinc fingers and is an important regulator of T-cell development and plays an important role in endothelial cell biology. Defects in this gene are the cause of hypoparathyroidism with sensorineural deafness and renal dysplasia. [provided by RefSeq, Nov 2009],disease:Defects in GATA3 are the cause of hypoparathyroidism with sensorineural deafness and renal dysplasia (HDR) [MIM:146255]; also known as Barakat syndrome.,function:Transcriptional activator which binds to the enhancer of the T-cell receptor alpha and delta genes. Binds to the consensus sequence 5'-AGATAG-3',similarity:Contains 2 GATA-type zinc fingers.,tissue specificity:T-cells and endothelial cells,


研究领域 (Research Area)

Immunology

图片 (Image Data)



Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using GATA3 (Phospho-Ser308) Antibody



Immunofluorescence analysis of HUVEC cells, using GATA3 (Phospho-Ser308) Antibody. The picture on the right is

**产品名称: GATA-3 (phospho Ser308) Rabbit Polyclonal
Antibody**
产品货号: APRab04726



blocked with the phospho peptide.

注意事项 (Note)

For research use only .