

产品名称: MAD1 (phospho Ser428) Rabbit Polyclonal Antibody
产品货号: APRab04966



产品概述 (Summary)

产品名称 (Production Name)	MAD1 (phospho Ser428) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAD1L1 MAD1L1; MAD1; TXBP181; Mitotic spindle assembly checkpoint protein
别名 (Alternative Names)	MAD1; Mitotic arrest deficient 1-like protein 1; MAD1-like protein 1; Mitotic checkpoint MAD1 protein homolog; HsMAD1; hMAD1; Tax-binding protein 181
基因 ID (Gene ID)	8379.0
蛋白 ID (SwissProt ID)	Q9Y6D9. The antiserum was produced against synthesized peptide derived from human MAD1 around the phosphorylation site of Ser428. AA range: 394-443

产品应用 (Application)

产品名称: MAD1 (phospho Ser428) Rabbit Polyclonal Antibody
产品货号: APRab04966



稀释比 (Dilution Ratio)

IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:5000-1:10000

蛋白分子量 (Molecular Weight)


研究背景 (Background)

MAD1L1 is a component of the mitotic spindle-assembly checkpoint that prevents the onset of anaphase until all chromosome are properly aligned at the metaphase plate. MAD1L1 functions as a homodimer and interacts with MAD2L1. MAD1L1 may play a role in cell cycle control and tumor suppression. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2015], disease: Defects in MAD1L1 are involved in the development and/or progression of various types of cancer., function: Component of the spindle-assembly checkpoint that prevents the onset of anaphase until all chromosomes are properly aligned at the metaphase plate. May recruit MAD2L1 to unattached kinetochores. Has a role in the correct positioning of the septum. Required for anchoring MAD2L1 to the nuclear periphery., induction: Increased by TP53., PTM: Phosphorylated; by BUB1. Become hyperphosphorylated in late S through M phases or after mitotic spindle damage. Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Belongs to the MAD1 family., subcellular location: From the beginning to the end of mitosis, it is seen to move from a diffusely nuclear distribution to the centrosome, to the spindle midzone and finally to the midbody., subunit: Homodimer. Heterodimerizes with MAD2L1 in order to form a tetrameric MAD1L1-MAD2L1 core complex. Perturbation of the original MAD1L1-MAD2L1 structure by the spindle checkpoint may decrease MAD2L1 affinity for MAD1L1. CDC20 can compete with MAD1L1 for MAD2L1 binding, until the attachment and/or tension dampen the checkpoint signal, preventing further release of MAD2L1 on to CDC20. Also able to interact with the BUB1/BUB3 complex and the viral Tax protein. Interacts with TPR., tissue specificity: Expressed weakly at G0/G1 and highly at late S and G2/M phase.,

研究领域 (Research Area)

Cell_Cycle_G1S; Cell_Cycle_G2M_DNA;

图片 (Image Data)




Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-

产品名称: MAD1 (phospho Ser428) Rabbit Polyclonal Antibody
产品货号: APRab04966



Phosphopeptide (Phospho-right) , using MAD1 (Phospho-Ser428) Antibody



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using MAD1 (Phospho-Ser428) Antibody.
The picture on the right is blocked with the phospho peptide.

注意事项 (Note)

For research use only .