

产品名称: Rrn3 (phospho Ser649) Rabbit Polyclonal Antibody
产品货号: APRab05390



产品概述 (Summary)

产品名称 (Production Name)	Rrn3 (phospho Ser649) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RRN3
别名 (Alternative Names)	RRN3; TIFIA; RNA polymerase I-specific transcription initiation factor RRN3; Transcription initiation factor IA; TIF-IA
基因 ID (Gene ID)	54700.0
蛋白 ID (SwissProt ID)	Q9NYV6. The antiserum was produced against synthesized peptide derived from human TIF-IA around the phosphorylation site of Ser649. AA range:602-651

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
----------------------	--

产品名称: Rrn3 (phospho Ser649) Rabbit Polyclonal Antibody
产品货号: AP Rab05390



蛋白分子量 (Molecular Weight) 70kDa


研究背景 (Background)

function: Required for efficient transcription initiation by RNA polymerase I., similarity: Belongs to the RRN3


family., function: Required for efficient transcription initiation by RNA polymerase I., similarity: Belongs to the RRN3 family.,

研究领域 (Research Area)

图片 (Image Data)




Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using TIF-IA (Phospho-Ser649) Antibody



Immunohistochemistry analysis of paraffin-embedded human brain, using TIF-IA (Phospho-Ser649) Antibody. The picture on the right is blocked with the phospho peptide.

产品名称: Rrn3 (phospho Ser649) Rabbit Polyclonal
Antibody
产品货号: APRab05390



Western blot analysis of lysates from Jurkat cells treated with starved 24h, using TIF-IA (Phospho-Ser649) Antibody. The lane on the right is blocked with the phospho peptide.

注意事项 (Note)

For research use only .