

产品名称: ZFP598 (phospho Tyr306) Rabbit Polyclonal Antibody
产品货号: APRab05646



产品概述 (Summary)

产品名称 (Production Name)	ZFP598 (phospho Tyr306) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZNF598
别名 (Alternative Names)	ZNF598; Zinc finger protein 598
基因 ID (Gene ID)	90850.0
蛋白 ID (SwissProt ID)	Q86UK7.Synthesized phospho-peptide around the phosphorylation site of human ZFP598 (phospho Tyr306)

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	99kDa

产品名称: ZFP598 (phospho Tyr306) Rabbit Polyclonal Antibody
产品货号: APRab05646




研究背景 (Background)


zinc finger protein 598(ZNF598) Homo sapiens Zinc-finger proteins bind nucleic acids and play important roles in various cellular functions, including cell proliferation, differentiation, and apoptosis. This protein and Grb10-interacting GYF protein 2 have been identified as a components of the mammalian 4EHP (m4EHP) complex. The complex is thought to function as a translation repressor in embryonic development. [provided by RefSeq, Oct 2012],similarity:Contains 1 C2H2-type zinc finger.,similarity:Contains 1 RING-type zinc finger.,

研究领域 (Research Area)

图片 (Image Data)



Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right) , using ZFP598 (Phospho-Tyr306) Antibody



Western blot analysis of lysate from Jurkat cells, using phospho-ZFP598 (Phospho-Tyr306) antibody.

注意事项 (Note)

For research use only .