

产品名称: PDE3A (Phospho-Ser312) Rabbit Polyclonal Antibody
产品货号: APRab05834



产品概述 (Summary)

产品名称 (Production Name)	PDE3A (Phospho-Ser312) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PDE3A
别名 (Alternative Names)	cGMP-inhibited 3',5'-cyclic phosphodiesterase A (EC 3.1.4.17) (Cyclic GMP-inhibited phosphodiesterase A) (CGI-PDE A)
基因 ID (Gene ID)	5139.0
蛋白 ID (SwissProt ID)	Q14432.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	126kDa

产品名称: PDE3A (Phospho-Ser312) Rabbit Polyclonal Antibody
产品货号: APRab05834



研究背景 (Background)

catalytic activity:Nucleoside 3',5'-cyclic phosphate + H(2)O = nucleoside 5'-phosphate.,enzyme regulation:Inhibited by cGMP.,function:Hydrolyzes both cyclic AMP (cAMP) and cyclic GMP (cGMP).,online information:PDE3 entry,similarity:Belongs to the cyclic nucleotide phosphodiesterase family.,

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of HeLa using primary antibody at 1:1000 dilution 4°C, overnight. Secondary antibody was diluted at 1:10000 25°C, 1.5hours

注意事项 (Note)

For research use only .