

产品名称: STAT5A/B (Phospho-Ser725/730) Rabbit Polyclonal Antibody
产品货号: APRab06088



产品概述 (Summary)

产品名称 (Production Name)	STAT5A/B (Phospho-Ser725/730) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	STAT5A
别名 (Alternative Names)	Signal transducer and activator of transcription 5A
基因 ID (Gene ID)	6776.0
蛋白 ID (SwissProt ID)	P42229.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	87kDa

研究背景 (Background)


产品名称: STAT5A/B (Phospho-Ser725/730) Rabbit
Polyclonal Antibody
产品货号: APRab06088



研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western Blot analysis of heLa 293T using primary antibody at 1:1000 dilution 4°C, overnight. Secondary antibody was diluted at 1:10000 25°C, 1.5hours

注意事项 (Note)

For research use only .