

产品名称: HNF-3 α / β / γ (Acetyl Lys264/253/211) Rabbit
Polyclonal Antibody
产品货号: APRab06218



产品概述 (Summary)

产品名称 (Production Name)	HNF-3 α / β / γ (Acetyl Lys264/253/211) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Acetylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FOXA1/FOXA2/FOXA3 FOXA1; HNF3A; TCF3A; Hepatocyte nuclear factor 3-alpha; HNF-3-alpha; HNF-3A; Forkhead box protein A1; Transcription factor 3A; TCF-3A; FOXA2; HNF3B; TCF3B; Hepatocyte nuclear factor 3-beta; HNF-3-beta; HNF-3B; Forkhead box protein A2; Transcription factor 3B; TCF-3B; FOXA3; HNF3G; TCF3G; Hepatocyte nuclear factor 3-gamma; HNF-3-gamma; HNF-3G; Fork head-related protein FKH H3; Forkhead box protein A3; Transcription factor 3G; TCF-3G
别名 (Alternative Names)	
基因 ID (Gene ID)	3169.0
蛋白 ID (SwissProt ID)	P55317. The antiserum was produced against synthesized Acetyl-peptide derived from human FOXA1 around the Acetylation site of Lys265. AA

产品名称: HNF-3 $\alpha/\beta/\gamma$ (Acetyl Lys264/253/211) Rabbit
Polyclonal Antibody
产品货号: APRab06218



range:231-280

产品应用 (Application)

稀释比 (Dilution Ratio) WB 1:500-1:2000, ELISA 1:10000-1:20000

蛋白分子量 (Molecular Weight) 50kDa


研究背景 (Background)

This gene encodes a member of the forkhead class of DNA-binding proteins. These hepatocyte nuclear factors are transcriptional activators for liver-specific transcripts such as albumin and transthyretin, and they also interact with chromatin. Similar family members in mice have roles in the regulation of metabolism and in the differentiation of the pancreas and liver. [provided by RefSeq, Jul 2008], function: Transcription activator for a number of liver genes such as AFP, albumin, tyrosine aminotransferase, PEPCCK, etc. Interacts with the cis-acting regulatory regions of these genes., online information: Hepatocyte nuclear factors entry, similarity: Contains 1 fork-head DNA-binding domain.,

研究领域 (Research Area)


Epigenetics and Nuclear Signaling

图片 (Image Data)




Western blot analysis of lysate from Jurkat cells, using FOXA1 (Acetyl-Lys265) Antibody.

产品名称: HNF-3 α / β / γ (Acetyl Lys264/253/211) Rabbit
Polyclonal Antibody
产品货号: APRab06218



Western Blot analysis of Jurkat cells using Acetyl-HNF-3 α / β / γ (K264/253/211) Polyclonal Antibody.. Secondary antibody was diluted at 1:20000



Western Blot analysis of JK cells using Acetyl-HNF-3 α / β / γ (K264/253/211) Polyclonal Antibody. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .