

产品名称: HSP 90 (Acetyl Lys435) Rabbit Polyclonal Antibody
产品货号: APRab06222



产品概述 (Summary)

产品名称 (Production Name)	HSP 90 (Acetyl Lys435) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Acetylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSP90AA1
别名 (Alternative Names)	HSP90AA1; HSP90A; HSPC1; HSPCA; Heat shock protein HSP 90-alpha; Heat shock 86 kDa; HSP 86; HSP86; Renal carcinoma antigen NY-REN-38
基因 ID (Gene ID)	3320.0
蛋白 ID (SwissProt ID)	P07900.Synthesized acetyl-peptide derived from the human HSP 90 around the acetylation site of K435.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	85kDa

产品名称: HSP 90 (Acetyl Lys435) Rabbit Polyclonal Antibody
产品货号: APRab06222




研究背景 (Background)

The protein encoded by this gene is an inducible molecular chaperone that functions as a homodimer. The encoded protein aids in the proper folding of specific target proteins by use of an ATPase activity that is modulated by co-chaperones. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2012],caution:Despite classification as a pseudogene, the existence of this protein is supported by unambiguous mass spectrometry evidence.,function:Molecular chaperone.,function:Molecular chaperone. Has ATPase activity.,similarity:Belongs to the heat shock protein 90 family.,subcellular location:Identified by mass spectrometry in melanosome fractions from stage I to stage IV.,subunit:Homodimer. Interacts with AHSA1, SMYD3 and TOM34. Interacts with FNIP1 and HSF1.,


研究领域 (Research Area)

Antigen processing and presentation;NOD-like receptor;Progesterone-mediated oocyte maturation;Pathways in cancer;Prostate cancer;

图片 (Image Data)



Western Blot analysis of A549 cells using Acetyl-HSP 90 (K435) Polyclonal Antibody.. Secondary antibody was diluted at 1:20000



Western Blot analysis of A549 cells using Acetyl-HSP 90 (K435) Polyclonal Antibody. Secondary antibody was diluted at 1:20000

产品名称: HSP 90 (Acetyl Lys435) Rabbit Polyclonal
Antibody
产品货号: APRab06222



注意事项 (Note)

For research use only .