

产品名称: 14-3-3 β/ζ Rabbit Polyclonal Antibody
产品货号: APRab06279



产品概述 (Summary)

产品名称 (Production Name)	14-3-3 β/ζ Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	YWHAB/YWHAZ YWHAZ; 14-3-3 protein zeta/delta; Protein kinase C inhibitor protein 1; KCIP-1;
别名 (Alternative Names)	YWHAB; 14-3-3 protein beta/alpha; Protein 1054; Protein kinase C inhibitor protein 1; KCIP-1
基因 ID (Gene ID)	7534/7529
蛋白 ID (SwissProt ID)	P63104/P31946.The antiserum was produced against synthesized peptide derived from human YWHAZ. AA range:151-200

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	


产品名称: 14-3-3 β/ζ Rabbit Polyclonal Antibody
产品货号: APRab06279



研究背景 (Background)

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from NIH/3T3 cells, using 14-3-3 beta/zeta Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .