

产品名称: 14-3-3 γ Rabbit Polyclonal Antibody
产品货号: APRab06280



产品概述 (Summary)

产品名称 (Production Name)	14-3-3 γ Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	YWHAG
别名 (Alternative Names)	YWHAG; 14-3-3 protein gamma; Protein kinase C inhibitor protein 1; KCIP-1
基因 ID (Gene ID)	7532.0
蛋白 ID (SwissProt ID)	P61981.The antiserum was produced against synthesized peptide derived from human 14-3-3 gamma. AA range:51-100

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)


产品名称: 14-3-3 γ Rabbit Polyclonal Antibody
产品货号: APRab06280




研究领域 (Research Area)

Neuroscience


图片 (Image Data)



Immunofluorescence analysis of COS7 cells, using 14-3-3 gamma Antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using 14-3-3 gamma Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from K562 cells, treated with insulin 0.01U/ml 15', using 14-3-3 gamma Antibody. The lane

产品名称: 14-3-3 γ Rabbit Polyclonal Antibody
产品货号: APRab06280



on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .