

**产品名称: AARSD1 Rabbit Polyclonal Antibody**  
**产品货号: APRab06380**



## 产品概述 (Summary)

|                        |                                   |
|------------------------|-----------------------------------|
| 产品名称 (Production Name) | AARSD1 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody        |
| 宿主 (Host)              | Rabbit                            |
| 应用 (Application)       | WB,ELISA                          |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                   |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                           |
|------------------------|-----------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | AARSD1                                                                                                    |
| 别名 (Alternative Names) | AARSD1; Alanyl-tRNA editing protein Aarsd1; Alanyl-tRNA synthetase domain-containing protein 1            |
| 基因 ID (Gene ID)        | 80755.0                                                                                                   |
| 蛋白 ID (SwissProt ID)   | Q9BTE6.The antiserum was produced against synthesized peptide derived from human AARSD1. AA range:141-190 |

## 产品应用 (Application)

|                          |                                       |
|--------------------------|---------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,ELISA 1:20000-1:40000 |
| 蛋白分子量 (Molecular Weight) | 45kDa                                 |

产品名称: AARSD1 Rabbit Polyclonal Antibody  
产品货号: APRab06380




## 研究背景 (Background)


caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,similarity:Belongs to the class-II aminoacyl-tRNA synthetase family. AARSD1 subfamily.,caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,similarity:Belongs to the class-II aminoacyl-tRNA synthetase family. AARSD1 subfamily.,

## 研究领域 (Research Area)

## 图片 (Image Data)




Western blot analysis of lysates from A549, 293, and HepG2 cells, using AARSD1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HUVEC cells using AARSD1 antibody.

产品名称: AARSD1 Rabbit Polyclonal Antibody  
产品货号: APRab06380



Western Blot analysis of various cells using AARSD1 Polyclonal Antibody

### 注意事项 (Note)

For research use only .