

产品名称: ABHD11 Rabbit Polyclonal Antibody
产品货号: APRab06432



产品概述 (Summary)

产品名称 (Production Name)	ABHD11 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ABHD11
别名 (Alternative Names)	ABHD11; WBSCR21; PP1226; Abhydrolase domain-containing protein 11; Williams-Beuren syndrome chromosomal region 21 protein
基因 ID (Gene ID)	83451.0
蛋白 ID (SwissProt ID)	Q8NFV4.The antiserum was produced against synthesized peptide derived from human ABHD11. AA range:161-210

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	32kDa

产品名称: ABHD11 Rabbit Polyclonal Antibody
产品货号: AP Rab06432




研究背景 (Background)


This gene encodes a protein containing an alpha/beta hydrolase fold domain. This gene is deleted in Williams syndrome, a multisystem developmental disorder caused by the deletion of contiguous genes at 7q11.23. [provided by RefSeq, Mar 2016], caution: It is uncertain whether Met-1 or Met-10 is the initiator., sequence caution: Wrong choice of frame., similarity: Belongs to the AB hydrolase superfamily., tissue specificity: Ubiquitously expressed.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from COLO and HUVEC cells, using ABHD11 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using ABHD11 antibody.

产品名称: ABHD11 Rabbit Polyclonal Antibody
产品货号: APRab06432



Western Blot analysis of various cells using ABHD11 Polyclonal Antibody

注意事项 (Note)

For research use only .