

产品名称: ABHD14A Rabbit Polyclonal Antibody
产品货号: APRab06435



产品概述 (Summary)

产品名称 (Production Name)	ABHD14A Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ABHD14A
别名 (Alternative Names)	ABHD14A; Abhydrolase domain-containing protein 14A
基因 ID (Gene ID)	25864.0
蛋白 ID (SwissProt ID)	Q9BUJ0.The antiserum was produced against synthesized peptide derived from human ABHD14A. AA range:201-250

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	30kDa

研究背景 (Background)


产品名称: ABHD14A Rabbit Polyclonal Antibody
产品货号: AP Rab06435



caution:It is uncertain whether Met-1 or Met-30 is the initiator.,function:Possible role in granule neuron development.,similarity:Belongs to the AB hydrolase superfamily. ABHD14 family.,caution:It is uncertain whether Met-1 or Met-30 is the initiator.,function:Possible role in granule neuron development.,similarity:Belongs to the AB hydrolase superfamily. ABHD14 family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from RAW264.7 cells, using ABHD14A Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .