

产品名称: ABHD2 Rabbit Polyclonal Antibody
产品货号: APRab06437



产品概述 (Summary)

产品名称 (Production Name)	ABHD2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ABHD2
别名 (Alternative Names)	ABHD2; LABH2; Abhydrolase domain-containing protein 2; Lung alpha/beta hydrolase 2; Protein PHPS1-2
基因 ID (Gene ID)	11057.0
蛋白 ID (SwissProt ID)	P08910. The antiserum was produced against synthesized peptide derived from human ABHD2. AA range: 251-300

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	48kDa

产品名称: ABHD2 Rabbit Polyclonal Antibody
产品货号: APRab06437




研究背景 (Background)


This gene encodes a protein containing an alpha/beta hydrolase fold, which is a catalytic domain found in a very wide range of enzymes. The function of this protein has not been determined. Alternative splicing of this gene results in two transcript variants encoding the same protein. [provided by RefSeq, Jul 2008], caution: Was originally (PubMed:2843827) thought to be a G-coupled receptor., function: May play a role in smooth muscle cells migration., similarity: Belongs to the AB hydrolase superfamily. AB hydrolase 4 family.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of MCF-7 cells, using ABHD2 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from HT29 cells, using ABHD2 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: ABHD2 Rabbit Polyclonal Antibody
产品货号: APRab06437



Western Blot analysis of various cells using ABHD2 Polyclonal Antibody

注意事项 (Note)

For research use only .