

产品名称: ABHD3 Rabbit Polyclonal Antibody
产品货号: APRab06438



产品概述 (Summary)

产品名称 (Production Name)	ABHD3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ABHD3 LABH3
别名 (Alternative Names)	
基因 ID (Gene ID)	171586.0
蛋白 ID (SwissProt ID)	Q8WU67.Synthesized peptide derived from human protein . at AA range: 240-320

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	44kDa

研究背景 (Background)

This gene encodes a protein containing an alpha/beta hydrolase fold, which is a catalytic domain found in a very wide


产品名称: ABHD3 Rabbit Polyclonal Antibody
产品货号: APRab06438



range of enzymes. The function of this protein has not been determined. [provided by RefSeq, Jul 2008],similarity:Belongs to the AB hydrolase superfamily. AB hydrolase 4 family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .