

产品名称: ABHD6 Rabbit Polyclonal Antibody
产品货号: APRab06441



产品概述 (Summary)

产品名称 (Production Name)	ABHD6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ABHD6
别名 (Alternative Names)	ABHD6; Monoacylglycerol lipase ABHD6; 2-arachidonoylglycerol hydrolase; Abhydrolase domain-containing protein 6
基因 ID (Gene ID)	57406.0
蛋白 ID (SwissProt ID)	Q9BV23. The antiserum was produced against synthesized peptide derived from human ABHD6. AA range: 231-280

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	38kDa

产品名称: ABHD6 Rabbit Polyclonal Antibody
产品货号: APRab06441




研究背景 (Background)


catalytic activity:Hydrolyzes glycerol monoesters of long-chain fatty acids.,function:Has 2-arachidonoylglycerol hydrolase activity (By similarity). May be a regulator of endocannabinoid signaling pathways.,similarity:Belongs to the AB hydrolase superfamily.,catalytic activity:Hydrolyzes glycerol monoesters of long-chain fatty acids.,function:Has 2-arachidonoylglycerol hydrolase activity (By similarity). May be a regulator of endocannabinoid signaling pathways.,similarity:Belongs to the AB hydrolase superfamily.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of A549 cells, using ABHD6 Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from COLO cells, using ABHD6 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: ABHD6 Rabbit Polyclonal Antibody
产品货号: APRab06441



Western Blot analysis of various cells using ABHD6 Polyclonal Antibody

注意事项 (Note)

For research use only .