

产品名称: ABHD9 Rabbit Polyclonal Antibody
产品货号: APRab06444



产品概述 (Summary)

产品名称 (Production Name)	ABHD9 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EPHX3
别名 (Alternative Names)	EPHX3; ABHD9; Epoxide hydrolase 3; Abhydrolase domain-containing protein 9
基因 ID (Gene ID)	79852.0
蛋白 ID (SwissProt ID)	Q9H6B9.The antiserum was produced against synthesized peptide derived from human ABHD9. AA range:101-150

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	38kDa

产品名称: ABHD9 Rabbit Polyclonal Antibody
产品货号: APRab06444




研究背景 (Background)


EPHX3 (Epoxide Hydrolase 3) is a Protein Coding gene. GO annotations related to this gene include hydrolase activity. An important paralog of this gene is EPHX4.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from MCF-7 cells, using ABHD9 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemical analysis of paraffin-embedded Human lung cancer. Antibody was diluted at 1:100 (4°, overnight). High-pressure and temperature Tris-EDTA, pH 8.0 was used for antigen retrieval. Negative control (right) obtained from antibody was pre-absorbed by immunogen peptide.

注意事项 (Note)

For research use only.