

产品名称: Abin-2 Rabbit Polyclonal Antibody
产品货号: APRab06450



产品概述 (Summary)

产品名称 (Production Name)	Abin-2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TNIP2
别名 (Alternative Names)	TNFAIP3-interacting protein 2 (A20-binding inhibitor of NF-kappa-B activation 2; ABIN-2; Fetal liver LKB1-interacting protein)
基因 ID (Gene ID)	79155.0
蛋白 ID (SwissProt ID)	Q8NFZ5. Synthetic peptide from human protein at AA range: 61-110

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:200, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)


产品名称: Abin-2 Rabbit Polyclonal Antibody
产品货号: APRab06450




This gene encodes a protein which acts as an inhibitor of NFkappaB activation. The encoded protein is also involved in MAP/ERK signaling pathway in specific cell types. It may be involved in apoptosis of endothelial cells. Alternative splicing results in multiple transcript variants. A pseudogene related to this gene is located on the X chromosome.[provided by RefSeq, May 2014],function:Inhibits NF-kappa-B activation by blocking the interaction of RIPK1 with its downstream effector IKBKG. Forms a ternary complex with NFKB1 and MAP3K8 but appears to function upstream of MAP3K8 in the TLR4 signaling pathway that regulates MAP3K8 activation.,subunit:Interacts with STK11, TNFAIP3, IKBKG, NFKB1, MAP3K8 and RIPK1.,tissue specificity:Ubiquitously expressed in all tissues examined.,

研究领域 (Research Area)

图片 (Image Data)




Immunohistochemical analysis of paraffin-embedded human-cervical-cancer, antibody was diluted at 1:200




Immunohistochemical analysis of paraffin-embedded human-cervical-cancer, antibody was diluted at 1:200

产品名称: Abin-2 Rabbit Polyclonal Antibody
产品货号: APRab06450



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200

注意事项 (Note)

For research use only .