

**产品名称: ADAM32 Rabbit Polyclonal Antibody**  
**产品货号: APRab06591**



## 产品概述 (Summary)

|                        |                                   |
|------------------------|-----------------------------------|
| 产品名称 (Production Name) | ADAM32 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody        |
| 宿主 (Host)              | Rabbit                            |
| 应用 (Application)       | WB,ELISA                          |
| 种属反应性 (Reactivity)     | Human,Rat,Mouse                   |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                    |
|------------------------|------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | ADAM32                                                                             |
| 别名 (Alternative Names) | ADAM32; Disintegrin and metalloproteinase domain-containing protein 32;<br>ADAM 32 |
| 基因 ID (Gene ID)        | 203102.0                                                                           |
| 蛋白 ID (SwissProt ID)   | Q8TC27.Synthesized peptide derived from ADAM32 . at AA range: 460-540              |

## 产品应用 (Application)

|                          |                                       |
|--------------------------|---------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,ELISA 1:20000-1:40000 |
| 蛋白分子量 (Molecular Weight) | 87kDa                                 |

## 研究背景 (Background)


**产品名称: ADAM32 Rabbit Polyclonal Antibody**  
**产品货号: APRab06591**



ADAM metallopeptidase domain 32(ADAM32) Homo sapiens This gene encodes a member of the disintegrin family of membrane-anchored proteins that play a role in diverse biological processes such as brain development, fertilization, tumor development and inflammation. This gene is predominantly expressed in the testis. The encoded protein undergoes proteolytic processing to generate a mature polypeptide comprised of an metalloprotease, disintegrin and epidermal growth factor-like domains. This gene is located in a cluster of other disintegrin and metallopeptidase family genes on chromosome 8. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2015],function:May play a role in sperm development and fertilization This is a non-catalytic metalloprotease-like protein.,similarity:Contains 1 disintegrin domain.,similarity:Contains 1 EGF-like domain.,similarity:Contains 1 peptidase M12B domain.,tissue specificity:Testis specific.,

## 研究领域 (Research Area)

## 图片 (Image Data)



Western Blot analysis of various cells using ADAM32 Polyclonal Antibody diluted at 1 : 1000

## 注意事项 (Note)

For research use only .