

产品名称: ADCK1 Rabbit Polyclonal Antibody
产品货号: APRab06608



产品概述 (Summary)

产品名称 (Production Name)	ADCK1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ADCK1
别名 (Alternative Names)	ADCK1; Uncharacterized aarF domain-containing protein kinase 1
基因 ID (Gene ID)	57143.0
蛋白 ID (SwissProt ID)	Q86TW2.The antiserum was produced against synthesized peptide derived from human ADCK1. AA range:251-300

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	61kDa

研究背景 (Background)


产品名称: ADCK1 Rabbit Polyclonal Antibody
产品货号: APRab06608




function: The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr)., similarity: Belongs to the protein kinase superfamily. ADCK protein kinase family., similarity: Contains 1 protein kinase domain., function: The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr)., similarity: Belongs to the protein kinase superfamily. ADCK protein kinase family., similarity: Contains 1 protein kinase domain.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from A549 and RAW264.7 cells, using ADCK1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using ADCK1 Polyclonal Antibody

注意事项 (Note)

For research use only .