

产品名称: AGR3 Rabbit Polyclonal Antibody
产品货号: AP Rab06680



产品概述 (Summary)

产品名称 (Production Name)	AGR3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AGR3
别名 (Alternative Names)	AGR3; BCMP11; Anterior gradient protein 3 homolog; AG-3; AG3; hAG-3; Breast cancer membrane protein 11
基因 ID (Gene ID)	155465.0
蛋白 ID (SwissProt ID)	Q8TD06. The antiserum was produced against synthesized peptide derived from human AGR3. AA range: 1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:200-1:1000, ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	

产品名称: AGR3 Rabbit Polyclonal Antibody
产品货号: AP Rab06680




研究背景 (Background)


similarity:Belongs to the AGR family.,subunit:Interacts with LYPD3 and alpha-dystroglycan.,similarity:Belongs to the AGR family.,subunit:Interacts with LYPD3 and alpha-dystroglycan.,

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of COS7 cells, using AGR3 Antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using AGR3 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .