

产品名称: AHNK Rabbit Polyclonal Antibody
产品货号: APRab06692



产品概述 (Summary)

产品名称 (Production Name)	AHNK Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC/IF
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AHNAK
别名 (Alternative Names)	PM227
基因 ID (Gene ID)	79026.0
蛋白 ID (SwissProt ID)	Q09666.Synthesized peptide derived from human protein . at AA range: 5720-5800

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:300, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	647kDa

研究背景 (Background)

function: May be required for neuronal cell differentiation., similarity: Contains 1 PDZ (DHR) domain., subunit: Interacts with


产品名称: AHNK Rabbit Polyclonal Antibody
产品货号: APRab06692



DYSF; the interaction is direct and Ca(2+)-independent.,function:May be required for neuronal cell differentiation.,similarity:Contains 1 PDZ (DHR) domain.,subunit:Interacts with DYSF; the interaction is direct and Ca(2+)-independent.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Tris-EDTA,pH9.0 was used for antigen retrieval. 2 Antibody was diluted at 1:200 (4° overnight.3,Secondary antibody was diluted at 1:200 (room temperature, 45min) .

注意事项 (Note)

For research use only .