

产品名称: AIG-1 Rabbit Polyclonal Antibody
产品货号: APRab06702



产品概述 (Summary)

产品名称 (Production Name)	AIG-1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AIG1
别名 (Alternative Names)	AIG1; CGI-103; Androgen-induced gene 1 protein; AIG-1
基因 ID (Gene ID)	51390.0
蛋白 ID (SwissProt ID)	Q9NVV5.The antiserum was produced against synthesized peptide derived from human AIG1. AA range:196-245

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	32kDa

研究背景 (Background)


产品名称: AIG-1 Rabbit Polyclonal Antibody
产品货号: APRab06702




function: May play a role in androgen-regulated growth of hair follicles., induction: By dihydrotestosterone (DHT)., similarity: Belongs to the AIG1 family., tissue specificity: Highly expressed in heart, ovary, testis, liver, and kidney, at lower levels in spleen, prostate, brain, skeletal muscle, pancreas, small intestine and colon, and undetected in peripheral blood leukocytes, thymus, lung and placenta. AIG1 expression is higher in hair follicles from males than from females., function: May play a role in androgen-regulated growth of hair follicles., induction: By dihydrotestosterone (DHT)., similarity: Belongs to the AIG1 family., tissue specificity: Highly expressed in heart, ovary, testis, liver, and kidney, at lower levels in spleen, prostate, brain, skeletal muscle, pancreas, small intestine and colon, and undetected in peripheral blood leukocytes, thymus, lung and placenta. AIG1 expression is higher in hair follicles from males than from females.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from HeLa and Jurkat cells, using AIG1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HT-29 cells using AIG1 antibody.

产品名称: AIG-1 Rabbit Polyclonal Antibody
产品货号: APRab06702



Western Blot analysis of various cells using AIG-1 Polyclonal Antibody

注意事项 (Note)

For research use only .