

产品名称: AKAP9 Rabbit Polyclonal Antibody
产品货号: APRab06730



产品概述 (Summary)

产品名称 (Production Name)	AKAP9 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC/IF
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AKAP9
别名 (Alternative Names)	AKAP350 AKAP450 KIAA0803
基因 ID (Gene ID)	10142.0
蛋白 ID (SwissProt ID)	Q99996.Synthesized peptide derived from human protein . at AA range: 2610-2690

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:300, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	430kDa

研究背景 (Background)

The A-kinase anchor proteins (AKAPs) are a group of structurally diverse proteins which have the common function of


产品名称: AKAP9 Rabbit Polyclonal Antibody
产品货号: APRab06730



binding to the regulatory subunit of protein kinase A (PKA) and confining the holoenzyme to discrete locations within the cell. This gene encodes a member of the AKAP family. Alternate splicing of this gene results in at least two isoforms that localize to the centrosome and the Golgi apparatus, and interact with numerous signaling proteins from multiple signal transduction pathways. These signaling proteins include type II protein kinase A, serine/threonine kinase protein kinase N, protein phosphatase 1, protein phosphatase 2a, protein kinase C-epsilon and phosphodiesterase 4D3. [provided by RefSeq, Aug 2008],disease:Defects in AKAP9 are the cause of long QT syndrome type 11 (LQT11) [MIM:611820]. Long QT syndromes are heart disorders characterized by a prolonged QT interval on the ECG and polymorphic ventricular arrhythmias. They cause syncope and sudden death in response to exercise or emotional stress. They can present with a sentinel event of sudden cardiac death in infancy.,domain:RII-binding site, predicted to form an amphipathic helix, could participate in protein-protein interactions with a complementary surface on the R-subunit dimer.,function:Binds to type II regulatory subunits of protein kinase A. Scaffolding protein that assembles several protein kinases and phosphatases on the centrosome and Golgi apparatus. May be required to maintain the integrity of the Golgi apparatus. Isoform 4/Yotiao is associated with the N-methyl-D-aspartate receptor and is specifically found in the neuromuscular junction (NMJ) as well as in neuronal synapses, suggesting a role in the organization of postsynaptic specializations.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the protein kinase superfamily.,subcellular location:Cytoplasmic in parietal cells.,subunit:Interacts with the regulatory region of protein kinase N (PKN), protein phosphatase 2A (PP2A), protein phosphatase 1 (PP1) and the immature non-phosphorylated form of PKC epsilon. Interacts with CIP4 and FNBP1.,tissue specificity:Widely expressed. Isoform 4/Yotiao is highly expressed in skeletal muscle and in pancreas.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight) . 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 30min) .

注意事项 (Note)

For research use only .

产品名称: AKAP9 Rabbit Polyclonal Antibody
产品货号: APRab06730

