

产品名称: ALKB7 Rabbit Polyclonal Antibody
产品货号: APRab06785



产品概述 (Summary)

产品名称 (Production Name)	ALKB7 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ALKBH7 ABH7 SPATA11 UNQ6002/PRO34564
别名 (Alternative Names)	
基因 ID (Gene ID)	84266.0
蛋白 ID (SwissProt ID)	Q9BT30.Synthesized peptide derived from human protein . at AA range: 10-90

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	24kDa

研究背景 (Background)


similarity:Belongs to the alkB family.,similarity:Belongs to the alkB family.,

产品名称: ALKB7 Rabbit Polyclonal Antibody
产品货号: APRab06785



研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only .