

产品名称: ALS Rabbit Polyclonal Antibody
产品货号: APRab06809



产品概述 (Summary)

产品名称 (Production Name)	ALS Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IGFALS
别名 (Alternative Names)	IGFALS; ALS; Insulin-like growth factor-binding protein complex acid labile subunit; ALS
基因 ID (Gene ID)	3483.0
蛋白 ID (SwissProt ID)	P35858.The antiserum was produced against synthesized peptide derived from human IGFALS. AA range:303-352

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	66kDa

产品名称: ALS Rabbit Polyclonal Antibody
产品货号: APRab06809




研究背景 (Background)


The protein encoded by this gene is a serum protein that binds insulin-like growth factors, increasing their half-life and their vascular localization. Production of the encoded protein, which contains twenty leucine-rich repeats, is stimulated by growth hormone. Defects in this gene are a cause of acid-labile subunit deficiency, which manifests itself in a delayed and slow puberty. Three transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Mar 2009],function:Involved in protein-protein interactions that result in protein complexes, receptor-ligand binding or cell adhesion.,similarity:Contains 20 LRR (leucine-rich) repeats.,subunit:Forms a ternary complex of about 140 to 150 kDa with IGF-I or IGF-II and IGFBP-3.,tissue specificity:Plasma.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of IGFALS Antibody. The lane on the right is blocked with the IGFALS peptide.



Western blot analysis of the lysates from RAW264.7 cells using IGFALS antibody.

产品名称: ALS Rabbit Polyclonal Antibody
产品货号: APRab06809



Western Blot analysis of CoLo cells using ALS Polyclonal Antibody

注意事项 (Note)

For research use only .